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



1. Mizushima N. (2007). Autophagy: process and function. *Genes Dev.* 21:2861–2873. doi: 10.1101/gad.1599207
2. Glick, D., Barth, S., & Macleod, K. F. (2010). Autophagy: cellular and molecular mechanisms. *The Journal of Pathology.* 221(1), 3–12. <http://doi.org/10.1002/path.2697>
3. Yamada, E., & Singh, R. (2012). Mapping Autophagy on to Your Metabolic Radar. *Diabetes.* 61(2), 272–280. <http://doi.org/10.2337/db11-1199>
4. van der Vaart, A., Mari, M., & Reggiori, F. (2008). A picky eater: exploring the mechanisms of selective autophagy in human pathologies. *Traffic.* 9(3), 281–9. doi: 10.1111/j.1600-0854.2007.00674.x
5. Mancias, J. D., & Kimmelman, A. C. (2016). Mechanisms of selective autophagy in normal physiology and cancer. *Journal of Molecular Biology,* 428(9 Pt A), 1659–1680. <http://doi.org/10.1016/j.jmb.2016.02.027>
6. Lippai M., & Lów P. (2014). “The Role of the Selective Adaptor p62 and Ubiquitin-Like Proteins in Autophagy,” *BioMed Research International.* 2014, Article ID 832704, 11 pages. doi:10.1155/2014/832704
7. Klionsky, D. J., Abdelmohsen, K., Abe, A., Abedin, et al. (2016). Guidelines for the use and interpretation of assays for monitoring autophagy (3rd edition). *Autophagy.* 12(1), 1–222. <http://doi.org/10.1080/15548627.2015.1100356>
8. Mizushima, N., Yoshimori, T., & Levine, B. (2010). Methods in Mammalian Autophagy Research. *Cell.* 140(3), 313–326. <http://doi.org/10.1016/j.cell.2010.01.028>
9. Vakifahmetoglu-Norberg, H., Xia, H., & Yuan, J. (2015). Pharmacologic agents targeting autophagy. *The Journal of Clinical Investigation,* 125(1), 5–13. <http://doi.org/10.1172/JCI73937>
10. Galluzzi L., Bravo-San Pedro J.M., Levine B., Green D.R., Kroemer G. (2017). Pharmacological modulation of autophagy: therapeutic potential and persisting obstacles. *Nat Rev Drug Discov.* 16(7), 487–511. <http://dx.doi.org/10.1038/nrd.2017.22>
11. Shoji-Kawata, S., Sumpter, R., Leveno, M., et al. (2013). Identification of a candidate therapeutic autophagy-inducing peptide. *Nature.* 494(7436), 201–206. <http://doi.org/10.1038/nature11866>
12. Füllgrabe J., Ghislat G., Cho D.H., Rubinsztein D.C. (2016). 'Transcriptional regulation of mammalian autophagy at a glance'. *Journal of Cell Science.* 129 (16), 3059–3066. doi: 10.1242/jcs.188920
13. Füllgrabe, J., Klionsky, D. J., & Joseph, B. (2014). The return of the nucleus: Transcriptional and epigenetic control of autophagy. *Nature Reviews Molecular Cell Biology.* 15(1), 65–74. doi: 10.1038/nrm3716
14. Lapierre, L. R., Kumsta, C., Sandri, M., Ballabio, A., & Hansen, M. (2015). Transcriptional and epigenetic regulation of autophagy in aging. *Autophagy.* 11(6), 867–880. <http://doi.org/10.1080/15548627.2015.1034410>
15. Füllgrabe, J., Lynch-Day, M. A., Heldring, N., et al. (2013). The histone H4 lysine 16 acetyltransferase hMOF regulates the outcome of autophagy. *Nature.* 500(7463), 468–471. doi: 10.1038/nature12313
16. Rohde, J. R., & Cardenas, M. E. (2003). The Tor Pathway Regulates Gene Expression by Linking Nutrient Sensing to Histone Acetylation. *Molecular and Cellular Biology.* 23(2), 629–635. <http://doi.org/10.1128/MCB.23.2.629-635.2003>
17. Khalil, H., Tazi, M., Caution, K., et al. (2016). Aging is associated with hypermethylation of autophagy genes in macrophages. *Epigenetics,* 11(5), 381–388. <http://doi.org/10.1080/15592294.2016.1144007>
18. Artal-Martinez de Narvajás, A., Gomez, T. S., Zhang, J.-S., Mann, A. O., et al. (2013). Epigenetic Regulation of Autophagy by the Methyltransferase G9a. *Molecular and Cellular Biology,* 33(20), 3983–3993. <http://doi.org/10.1128/MCB.00813-13>

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<b>AMBRA1</b>	Activating molecule BECN1 regulated autophagy protein 1	<b>NDP52</b>	Nuclear dot protein 52 kDa
<b>AMPK</b>	AMP-activated protein kinase	<b>NIX</b>	NIP3-like protein X
<b>ATG</b>	Autophagy related proteins	<b>OPTN</b>	Optineurin
<b>BME</b>	2-mercaptoethanol	<b>PB1</b>	Phox and Bem1 domain
<b>BNIP3</b>	BCL-2/adenovirus E1B 19 kDa interacting protein 3	<b>PE</b>	Phosphatidylethanolamine
<b>CD40</b>	Cluster of differentiation 40	<b>PKB</b>	Protein kinase B
<b>Cue5</b>	Coupling of ubiquitin conjugation to ER degradation	<b>PTD</b>	Protein transduction domain
<b>DNMT</b>	DNA methyltransferase	<b>PtdIns3K</b>	Class III phosphatidylinositol 3-kinase
<b>DOR</b>	Diabetes- and obesity-regulated	<b>RB1CC1</b>	RB1-inducible coiled-coil 1
<b>ELISA</b>	Enzyme linked immunosorbent assay	<b>RPD3</b>	Reduced Potassium Deficiency 3/Histone Deacetylase
<b>ESA</b>	Essential SAS2-related acetyltransferase	<b>RPN10</b>	Regulatory particle non-ATPase
<b>FAM134B</b>	Family with sequence similarity 134 member B	<b>SIRT1</b>	Silent Mating Type Information Regulation 2, <i>S. Cerevisiae</i> , Homolog 1
<b>FC</b>	Flow cytometry	<b>SQSTM1</b>	Sequestosome 1
<b>FIP200</b>	FAK family-interacting protein of 200 kDa	<b>SMURF1</b>	SMAD ubiquitylation regulatory factor
<b>FM</b>	Fluorescent microscopy	<b>SNAP29</b>	Synaptosome associated protein 29
<b>FOXO</b>	Forkhead box O	<b>STX17</b>	Syntaxin 17
<b>FUNDC1</b>	FUN-domain containing protein	<b>TAX1BP1</b>	Human T-Cell Leukemia Virus Type I) Binding Protein 1
<b>G9A</b>	Euchromatic histone-lysine N-methyltransferase 2	<b>TCR</b>	T-cell receptor
<b>GAPR-1</b>	Golgi-associated plant pathogenesis-related protein 1	<b>TEM</b>	Transmission electron microscopy
<b>HSC-70</b>	Heat shock cognate 70 protein	<b>TFEB</b>	Transcription factor EB
<b>Hsp</b>	Heat shock protein	<b>TLR4</b>	Toll-like receptor 4
<b>ICC</b>	Immunocytochemistry	<b>TOLLIP</b>	Toll-interacting protein
<b>IF</b>	Immunofluorescence	<b>TRAIL</b>	Tumor necrosis factor-related apoptosis-inducing ligand
<b>IHC</b>	Immunohistochemistry	<b>TRIM5a</b>	Tripartite motif
<b>ISH</b>	in situ hybridization	<b>UBD</b>	Ubiquitin binding domain
<b>KO</b>	Knockout	<b>ULK1</b>	unc-51-like kinase 1
<b>LAMP-2A</b>	Lysosomal associated membrane protein-2A	<b>USP</b>	Ubiquitin Specific Peptidase
<b>LC3</b>	Microtubule-associated protein light chain 3	<b>UVRAG</b>	UV resistance-associated gene
<b>LIR</b>	LC3 interacting region	<b>VAMP8</b>	Vesicle associated membrane protein 8
<b>MHC</b>	Major histocompatibility complex	<b>VPS</b>	Vacuolar sorting protein
<b>mTOR</b>	Mammalian target of rapamycin	<b>WB</b>	Western Blot
<b>NBR1</b>	neighbor of BRCA1 gene 1	<b>WIPI</b>	WD-repeat protein interacting with phosphoinositides
		<b>ZKSCAN3</b>	zinc finger with KRAB and SCAN domains 3

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